

# PROFIT FROM the Award-winning Curing Light

TIME SAVING - \$26,399 PER YEAR  
Generate New Revenue Stream - Oral Exam!



### Dr. Mark S. Geller, DDS, MSD, P.A

"My staff and I have used the Fusion curing light for several years and are very pleased with its functional design and durability. Its short curing time, due to a high powered focused beam, ensures complete adhesive set in minimal time that adds to both efficiency and productivity. My office highly endorses this product."

### Order Information

7500080-08 FUSION 3.0 Silver System:  
Including solid-metal wand, tacking tip, 9mm and 13mm  
tip, Battery Assembly, Eyeshield, Power Adapter,  
Charging Base, 100 Barrier Sleeves  
2-year warranty (light unit)

### DENTLIGHT INC

1411 E. Campbell Rd. #500  
Richardson, TX 75081  
Tel: 800-763-6901  
Fax: 972-346-6550  
Email: sales@dentlight.com  
www.dentlight.com

**DentLight**

CALL 800-763-6901 today

# FINALLY, the LIGHT THAT YOU CAN COUNT ON!



## FUSION™ Curing Light

Lightning Power and Speed!

V 3.0



**DentLight**

## The Most Reliable and Versatile Curing Light

### Amazing Cordless Power in the Palm of Your Hand!

The award winning and clinically proven cordless FUSION 3.0 outputs a laser-sharp focused beam of 2700 mW/cm<sup>2</sup> power that enables rapid cures in 5 seconds. The new release improves on ease of use and work efficiency with automatic power load and three preset digital timers for single layer, bulk, multiple restorations, and veneer cures with a easy one-push activation. Its amazing power, incredible engineering design and proven quality will make this curing light your office favorite.

- Cordless with rechargeable smart lithium-ion battery
- Unique focused beam with two power level settings
  - Plasma-like High Power mode for fast/deep cure
  - Standard LED Power mode for low heat cure
- Wide spectrum of 415 - 490nm to cure all materials
- 45 minute run time
  - 540 5-second bondings
- Shock-proof durable solid metal construction
  - No fiber optic light guide necessary
  - Easy cleaning and disinfection
  - Modular with easy upgrade and service parts
- Raised metal keypad with LCD digital displays
- Cordless with direct power plug that eliminates downtime
- 90 minute super fast recharge time
- 360° swivel head




### Built to Last

Robust anodized aluminum metal construction distributes LED heat, prevent scratches, discoloration and damage from accidental falls.

### Incredible Engineering Design

- Slim and ergonomic
  - Weighs only 3.9oz
  - Measures at 7.8 inch
- Sleek metal keypad with easy activation raised buttons
- Three digital presets for rapid cure of any procedures
- Battery bar indicator
- Cordless, pocket-size portability
- **Beyond adhesive cure (optional)**
  - High contrast fluorescence imaging for composite removal and oral exam including oral cancer screening
  - Transillumination for caries and lesions detection



Violet Light

White Light

### Low Profile

360° swivel head for easier maneuvering and access



FUSION™

Demetron™ LED II

### Non-memory Li-Ion Battery You can Count On!


- Cordless power that outshines others
- Multiple charge options
  - Space saving counter stand
  - Direct plug in
  - Chair side charge holder
- Battery status always monitored on the LCD screen
- Secondary charge indicator
- Constant output throughout battery run time

**Attract New Patient**  
-> Add Oral Screening

## The Most Advanced Curing Light


### Patented Laser-sharp Focused Beam

The award winning FUSION outputs an unmatched power of 2700 mW/cm<sup>2</sup> with an unparalleled uniform beam profile and the lowest dispersion



Ref. Richard B. Price; Daniel Labrie; J. Marc Whalen; Christopher M. Felix, "Effect of Distance on Irradiance and Beam Homogeneity from 4 Light-Emitting Diode Curing Units", J. Can Dent Assoc 2011; 77; b9

### The Best Curing Light and the Best Bang for Your Buck!



### Incredible size and unmatched durability

- Solid metal construction with shock proof reliability
- Modular, upgradable functions and easy service
- Cordless or corded with replaceable battery

### You've got the POWER!

- Plasma-equivalent cure speed
- 10-sec. cures harder than other lights in 40 seconds

Data courtesy of Reality Publishing - The Ratings for composite hardness (Knoop) tested at 2mm using the highest power tip from each brand. FUSION is a registered brand of DentLight. Other brands are not registered brands of DentLight.



Tacking Tip enables best crown/veneer contact with pressure tack

\* Photo Courtesy of S. Benjamin on Cerec Crown

### Amazing LCD digital timer in a light this small!



- 5, 10 and 20 digital timer modes preset
- Two power levels for plasma speed or low heat LED bonding
- One click cure with no annoying sleep mode
- Remember the last setting